

Rabbit Anti-Human/Monkey Cd20 Monoclonal Antibody For Fc - P11836

Item Number: CM0012534



Introduction

Rabbit anti-Human/Monkey CD20 monoclonal antibody, reactive with human and cynomolgus. Validated for flow cytometry, targeting B-lymphocyte antigen CD20, essential for B-cell development, differentiation, and activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	B1, S7, Bp35, CD20, FMC7, CVID5, LEU-16
Species	Human
Gene ID	931 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 (NP_068769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag