

## Rabbit Anti-Human/Monkey Cd20 Monoclonal Antibody For Flow Cytometry - P11836

Item Number: CM0011070



### Introduction

High-quality Rabbit anti-Human/Monkey CD20 monoclonal antibody validated for flow cytometry. Detects CD20, a B-cell-specific calcium channel (P11836), with cross-reactivity to human and cynomolgus monkey.

[Learn More](#)

| Parameter        | Value                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Monkey CD20 mAb                                                                        |
| Remarks/Alias    | B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16                                                                  |
| Species          | Human                                                                                                    |
| Gene ID (Human)  | 931                                                                                                      |
| Gene ID          | 931                                                                                                      |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 mAb (NP_068769.2) |
| Source           | Rabbit                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                    |
| Application      | FC                                                                                                       |
| Cross Reactivity | Human, Cynomolgus                                                                                        |
| SWISS            | P11836                                                                                                   |
| Protein Weight   | 14kDa/33kDa                                                                                              |
| Shipping         | Ice bag                                                                                                  |