


Rabbit Anti-Human/Monkey Cd20 Monoclonal Antibody For Flow Cytometry - P11836

Item Number: CM0006126



Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human/monkey CD20 antigen. Validated for flow cytometry. B-cell marker for immunology research, diagnostics development, and drug screening. UniProt: P11836.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	B1; S7; Bp35; CD20; FMC7; CVID5; LEU-16
Species	Human
GeneID (Human)	931
GeneID	931
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 (NP_068769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag