

Rabbit Anti-Human/Monkey Cd25 Monoclonal Antibody For Fc - P01589

Item Number: CM0011312



Introduction

-conjugated rabbit monoclonal antibody targeting CD25 (IL2RA). Validated for flow cytometry (FC) in human and cynomolgus. Ideal for studying regulatory T cells (TREGs) and immune tolerance.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
Gene ID	3559
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag