

Rabbit Anti-Human/Monkey Cd25 Monoclonal Antibody For Flow Cytometry - P01589

Item Number: CM0009124



Introduction

-conjugated rabbit monoclonal antibody against human/monkey CD25 (IL-2 receptor alpha). Validated for flow cytometry, suitable for detecting regulatory T cells and studying immune tolerance mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
Gene ID (Human)	3559
Gene ID	3559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-213 of human CD25 (NP_000408.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag