

Rabbit Anti-Human/Monkey Cd3 Monoclonal Antibody For Fc - P07766

Item Number: CM0004791



Introduction

Rabbit anti-Human/Monkey CD3 Monoclonal Antibody for flow cytometry. Detects CD3 epsilon chain (23 kDa) on T cells. Cross-reacts with human and cynomolgus. Suitable for immunophenotyping and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag