


Rabbit Anti-Human/Monkey Cd38 Monoclonal Antibody For Fc - P28907

Item Number: CM0009079



Introduction

-conjugated rabbit anti-human/monkey CD38 monoclonal antibody for flow cytometry. Targets the 14-34 kDa transmembrane glycoprotein on lymphocytes and other tissues. Ideal for immunophenotyping and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (NP_001766.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	14kDa/34kDa
Shipping	Ice bag