


Rabbit Anti-Human/Monkey Cd38 Monoclonal Antibody For Flow Cytometry - P28907

Item Number: CM0012673



Introduction

-conjugated rabbit monoclonal antibody against human and cynomolgus CD38, validated for flow cytometry. Ideal for studying CD38's enzymatic and receptor functions in immune regulation, diabetes, and oncology research.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human/Monkey CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (NP_001766.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	14kDa/34kDa
Shipping	Ice bag