


Rabbit Anti-Human/Monkey Cd38 Polymab® For Fc - P28907

Item Number: CM0008247



Introduction

Rabbit anti-Human/Monkey CD38 PolymAb® is a monoclonal antibody for flow cytometry, validated for human and cynomolgus monkey. Targets CD38 glycoprotein involved in immune cell signaling and metabolism.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD38 PolymAb®
Remarks/Alias	CD38; ADPRC 1; ADPRC1; CD38 molecule
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag