

Rabbit Anti-Human/Monkey Cd4 Mab For Flow Cytometry - P01730

Item Number: CM0011326



Introduction

Rabbit anti-Human/Monkey CD4 monoclonal antibody (conjugate) validated for flow cytometry. Detects CD4, a key T-helper cell marker and HIV receptor. Suitable for immune cell phenotyping and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
Gene ID (Human)	920
Gene ID	920
Immunogen	Recombinant protein of human CD4.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag