

Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Flow Cytometry - P01730

Item Number: CM0004840



Introduction

Rabbit monoclonal antibody against human/monkey CD4, conjugated to , validated for flow cytometry (FC). Targets T-cell surface glycoprotein CD4; reacts with human and cynomolgus. SWISS: P01730.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID	920
Immunogen	Recombinant protein of human CD4.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag