

Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Flow Cytometry - P01730

Item Number: CM0009939



Introduction

Rabbit anti-Human/Monkey CD4 mAb, validated for flow cytometry (FC). Targets human CD4 (P01730), cross-reacts with cynomolgus. Useful for T-helper cell identification and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag