

Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Flow Cytometry - P01730

Item Number: CM0005368



Introduction

This rabbit monoclonal antibody specifically recognizes human, cynomolgus, and rhesus CD4 protein. Conjugated with , it is validated for flow cytometry, suitable for T-helper cell identification and HIV research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
Gene ID (Human)	920
Gene ID	920
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus, Rhesus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag