

Rabbit Anti-Human/Monkey Cd45 Mab For Flow Cytometry - P08575-9

Item Number: CM0011914

Introduction

-conjugated rabbit anti-human/monkey CD45 monoclonal antibody for flow cytometry. Targets leukocyte common antigen, cross-reacts with human and cynomolgus. Suitable for immunophenotyping. Isoform 7, ~138 kDa.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575-9
Protein Weight	138kDa
Shipping	Ice bag