

Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Fc - P08575

Item Number: CM0004608



Introduction

-conjugated rabbit anti-human/monkey CD45 monoclonal antibody for flow cytometry. Recognizes the common leukocyte antigen (LCA). Suitable for immune cell profiling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	138kDa
Shipping	Ice bag