

Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575

Item Number: CM0008762



Introduction

-conjugated rabbit anti-CD45 monoclonal antibody for flow cytometry, cross-reactive with human and cynomolgus CD45. Detects leukocyte common antigen for immunophenotyping and T-cell activation studies.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA, LY5, B220, CD45, L-CA, T200, CD45R, GP180, IMD105
Species	Human
Gene ID	5788 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	~147 kDa (canonical; multiple isoforms detected at 131-147 kDa)
Shipping	Ice bag