

Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0010932

Introduction

Recombinant rabbit anti-CD45 monoclonal antibody conjugated to for flow cytometry. Reacts with human and cynomolgus monkey CD45 isoform 6 (P08575-8), ideal for immune cell profiling and phenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag