

Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575-9

Item Number: CM0009933

Introduction

-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD45 (PTPRC, LCA). Designed for flow cytometry (FC), this antibody recognizes isoform 7 (P08575-9, ~138kDa). Supplied in an ice-packed format for stability.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575-9
Protein Weight	138kDa
Shipping	Ice bag