

Rabbit Anti-Human/Monkey Cd45 Monoclonal Antibody For Flow Cytometry - P08575

Item Number: CM0005142



Introduction

-conjugated rabbit monoclonal antibody targeting human and monkey CD45, validated for flow cytometry. Useful for leukocyte common antigen detection in immunology research.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575
Protein Weight	147kDa/143kDa/141kDa/136kDa/131kDa
Shipping	Ice bag