

Rabbit Anti-Human/Monkey Cd64 Mab For Flow Cytometry - P12314

Item Number: CM0012189



Introduction

-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD64 (FcγRI). Verified for flow cytometry applications, suitable for monocyte/macrophage identification and immune function assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD64 mAb
Remarks/Alias	CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcγRI
Species	Human
Gene ID (Human)	2209
Gene ID	2209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-288 of human CD64 (NP_000557.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P12314
Protein Weight	43kDa
Shipping	Ice bag