

Rabbit Anti-Human/Monkey Cd64 Monoclonal Antibody For Flow Cytometry - P12314

Item Number: CM0009850



Introduction

-conjugated rabbit monoclonal antibody targeting human/monkey CD64 (FcγRI), validated for flow cytometry. Suitable for innate and adaptive immunity research. Cross-reactivity with cynomolgus. Monocyte/macrophage specific high affinity IgG receptor.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD64 mAb
Remarks/Alias	CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcgammaRI
Species	Human
GeneID (Human)	2209
GeneID	2209
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 16-288 of human CD64 (NP_000557.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P12314
Protein Weight	43kDa
Shipping	Ice bag