


Rabbit Anti-Human/Monkey Cd69 Monoclonal Antibody For Flow Cytometry - Q07108

Item Number: CM0007808



Introduction

-conjugated rabbit monoclonal antibody against human CD69 (SWISS: Q07108). Validated for flow cytometry on human and cynomolgus samples. Suitable for immune cell activation and homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID	969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag