


Rabbit Anti-Human/Monkey Cd86 Monoclonal Antibody For Fc - P42081

Item Number: CM0008050



Introduction

-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD86, validated for flow cytometry. CD86 is a costimulatory molecule on antigen-presenting cells essential for T-cell activation. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD86 mAb
Remarks/Alias	B70; B7-2; B7.2; LAB72; CD28LG2
Species	Human
Gene ID	942 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-241 of human CD86 (NP_787058.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P42081
Protein Weight	38kDa
Shipping	Ice bag