


## Rabbit Anti-Human/Monkey Cd8A Monoclonal Antibody For Fc - P01732

Item Number: CM0012388



### Introduction

-conjugated rabbit monoclonal antibody against human and cynomolgus CD8a, suitable for flow cytometry. Recognizes the T-cell coreceptor involved in MHC class I-mediated immune responses. SWISS: P01732.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Monkey CD8a mAb                                                                                                     |
| Remarks/Alias    | CD8; p32; Leu2; CD8alpha                                                                                                              |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 925                                                                                                                                   |
| Gene ID          | 925                                                                                                                                   |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8A (NP_001759.3) |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC                                                                                                                                    |
| Cross Reactivity | Human, Cynomolgus                                                                                                                     |
| SWISS            | P01732                                                                                                                                |
| Protein Weight   | 21kDa/25kDa/30kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                               |