

Rabbit Anti-Human/Monkey Hla-Dr Mab For Fc - P01903

Item Number: CM0011462



Introduction

-conjugated rabbit monoclonal antibody against HLA-DRA (MHC II DR alpha chain) for flow cytometry analysis of human and cynomolgus monkey immune cells. Ideal for immunophenotyping and immune response research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
GeneID (Human)	3122
GeneID	3122
Immunogen	Recombinant protein of human HLA-DRA.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag