

Rabbit Anti-Human/Monkey Hla-Dr Monoclonal Antibody For Flow Cytometry - P01903

Item Number: CM0012592



Introduction

-conjugated rabbit monoclonal antibody against human/monkey HLA-DR alpha chain. Validated for flow cytometry; ideal for studying MHC class II antigen presentation in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
Gene ID (Human)	3122
Gene ID	3122
Immunogen	Recombinant protein of human HLA-DRA
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag