

Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Fc - P01880

Item Number: CM0005628



Introduction

-conjugated rabbit anti-human/monkey IgD monoclonal antibody for flow cytometry (FC). Validated in human and rhesus. Suitable for B cell immunophenotyping and studying IgD-mediated immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgD mAb
Species	Human
GeneID (human)	3495
GeneID	3495
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Rhesus
SWISS	P01880
Protein Weight	42kDa
Shipping	Ice bag