


## Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Flow Cytometry - P01880

Item Number: CM0012776



### Introduction

Rabbit anti-Human/Monkey IgD monoclonal antibody validated for flow cytometry. Recognizes human and cynomolgus IgD constant region. Ideal for B-cell receptor studies and humoral immunity research. 47 kDa.

[Learn More](#)

| Parameter        | Value                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Monkey IgD mAb                                                                  |
| Remarks/Alias    | IGHD                                                                                              |
| Species          | Human                                                                                             |
| GeneID (Human)   | 3495                                                                                              |
| GeneID           | 3495                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD |
| Source           | Rabbit                                                                                            |
| Category         | Monoclonal Antibodies                                                                             |
| Application      | FC                                                                                                |
| Cross Reactivity | Human, Cynomolgus                                                                                 |
| SWISS            | P01880                                                                                            |
| Protein Weight   | 47kDa                                                                                             |
| Shipping         | Ice bag                                                                                           |