

Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Flow Cytometry - P01880

Item Number: CM0009656



Introduction

Rabbit anti-Human/Monkey IgD monoclonal antibody for Flow Cytometry. Recognizes human and cynomolgus IgD. Suitable for B-cell receptor studies, immunoglobulin research, and immunology applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgD mAb
Remarks/Alias	IGHD
Species	Human
Gene ID	3495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD.
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human, Cynomolgus
SWISS	P01880
Protein Weight	47kDa
Shipping	Ice bag