


Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Flow Cytometry - P01880

Item Number: CM0006250



Introduction

-conjugated rabbit monoclonal antibody targeting human/monkey immunoglobulin heavy constant delta (IgD) for flow cytometry. Validated for human and rhesus samples, with a protein weight of 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgD mAb
Species	Human
GeneID (Human)	3495
GeneID	3495
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Rhesus
SWISS	P01880
Protein Weight	42kDa
Shipping	Ice bag