

Rabbit Anti-Human/Monkey IgG Mab - Aaq57561.1

Item Number: CM0008768

Introduction

-conjugated rabbit monoclonal antibody targeting human and monkey IgG heavy chain. Designed for intracellular flow cytometry (FC intra) applications. Shows cross-reactivity with human and cynomolgus monkey IgG.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human/Monkey IgG mAb
Species	Monkey
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-322 of monkey IgG (AAQ57561.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag