

Rabbit Anti-Human/Monkey Igm Mab - P01871-1

Item Number: CM0007850

Introduction

-conjugated rabbit monoclonal antibody targeting the human/monkey IgM heavy chain constant region (P01871-1). Validated for IF/ICC and flow cytometry. Ideal for B-cell and immunoglobulin research in diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	IGHM; AGM1; MU; VH; immunoglobulin heavy constant mu
Species	Human
GeneID (Human)	3507
GeneID	3507
Immunogen	Recombinant fusion protein containing amino acids 218-453 of human IgM (P01871)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871-1
Protein Weight	49kDa
Shipping	Ice bag