

Rabbit Anti-Human/Monkey IgM Mab - P01871-1

Item Number: CM0004538

Introduction

-conjugated recombinant rabbit monoclonal antibody targeting human and cynomolgus IgM heavy chain constant region. Suitable for immunofluorescence (IF/ICC, IF-P) and flow cytometry applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	IGHM; AGM1; MU; VH; immunoglobulin heavy constant mu
Species	Human
Gene ID (Human)	3507
Gene ID	3507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-453 of human IgM (P01871)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871-1
Protein Weight	49kDa/51kDa
Shipping	Ice bag