


Rabbit Anti-Human/Monkey IgM Monoclonal Antibody For Fc - P01871

Item Number: CM0009051



Introduction

-conjugated rabbit monoclonal antibody specific for the constant region of human IgM, validated for flow cytometry (FC). Cross-reacts with cynomolgus monkey, ideal for B cell immunophenotyping and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	MU; VH; AGM1
Species	Human
Gene ID	3507
Immunogen	Recombinant fusion protein containing amino acids 218-453 of human IgM (P01871)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871
Protein Weight	49 kDa
Shipping	Ice bag