


Rabbit Anti-Human/Monkey IgM Monoclonal Antibody For Flow Cytometry - P01871

Item Number: CM0012281



Introduction

-conjugated rabbit monoclonal antibody targeting human/monkey IgM heavy chain constant region (P01871), designed for flow cytometry. Ideal for B cell and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	MU; VH; AGM1
Species	Human
GeneID (Human)	3507
GeneID	3507
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 218-453 of human IgM (P01871).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871
Protein Weight	49kDa
Shipping	Ice bag