

Rabbit Anti-Human/Monkey IgM Monoclonal Antibody For If/Icc & Fc - P01871-1

Item Number: CM0011097

Introduction

High-purity rabbit monoclonal antibody against human/cynomolgus IgM, conjugated to , validated for immunofluorescence and flow cytometry. Targets secreted IgM constant region (SWISS: P01871-1). Ideal for humoral immunity research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	IGHM; AGM1; MU; VH; immunoglobulin heavy constant mu
Species	Human
GeneID	3507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-453 of human IgM (P01871)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871-1
Protein Weight	49kDa
Shipping	Ice bag