


Rabbit Anti-Human/Monkey Ki67 Mab - P46013

Item Number: CM0012228



Introduction

Rabbit anti-Human/Monkey Ki67 monoclonal antibody designed for intracellular flow cytometry. Targets human Ki67 proliferation marker; cross-reacts with cynomolgus. Ideal for cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey Ki67 mAb
Remarks/Alias	KIA; MIB-; MIB-1; PPP1R105
Species	Human
GeneID	4288
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Cynomolgus
SWISS	P46013
Protein Weight	319kDa/359kDa
Shipping	Ice bag