

Rabbit Anti-Human/Monkey Ki67 Monoclonal Antibody For Flow Cytometry (Intracellular) - P46013

Item Number: CM0010317



Introduction

-conjugated rabbit monoclonal antibody against human and monkey Ki67 (P46013). Designed for intracellular flow cytometry, enabling direct detection of proliferating cells in research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey Ki67 mAb
Remarks/Alias	KIA; MIB-; MIB-1; PPP1R105
Species	Human
Gene ID (Human)	4288
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular) (FC (intra))
Cross Reactivity	Human, Cynomolgus
SWISS	P46013
Protein Weight	319kDa/359kDa
Shipping	Ice bag