

## Rabbit Anti-Human/Monkey Ki67 Monoclonal Antibody For Flow Cytometry - P46013

Item Number: CM0011365



### Introduction

Fluorophore-conjugated monoclonal antibody against human Ki67, validated for intracellular flow cytometry. Cross-reacts with cynomolgus monkey. Suitable for proliferation and cell cycle studies.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Monkey Ki67 mAb                                                                      |
| Remarks/Alias    | KIA; MIB-; MIB-1; PPP1R105                                                                             |
| Species          | Human                                                                                                  |
| Gene ID (Human)  | 4288                                                                                                   |
| Gene ID          | 4288                                                                                                   |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3) |
| Source           | Rabbit                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                  |
| Application      | FC (intra)                                                                                             |
| Cross Reactivity | Human, Cynomolgus                                                                                      |
| SWISS            | P46013                                                                                                 |
| Protein Weight   | 319kDa/359kDa                                                                                          |
| Shipping         | Ice bag                                                                                                |