

Rabbit Anti-Human/Monkey Ki67 Monoclonal Antibody For Intracellular Flow Cytometry - P46013

Item Number: CM0004283

Introduction

-conjugated rabbit monoclonal antibody targeting Ki67, a proliferation marker. Validated for intracellular flow cytometry in human and cynomolgus samples. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey Ki67 mAb
Remarks/Alias	KIA; MIB-; MIB-1; PPP1R105
Species	Human
GeneID (Human)	4288
GeneID	4288
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Cynomolgus
SWISS	P46013
Protein Weight	319kDa/359kDa
Shipping	Ice bag