

## Rabbit Anti-Human/Mouse Cd276/B7-H3 Monoclonal Antibody For Fc - Q5Zpr3-2

Item Number: CM0011499

### Introduction

-conjugated rabbit anti-human/mouse CD276 (B7-H3) monoclonal antibody validated for flow cytometry. Recognizes isoform 2 of the immune costimulatory molecule; predicted MW ~34 kDa, observed at 33/52/57 kDa. Cross-reacts with human and mouse.

[Learn More](#)

| Parameter        | Value                                                                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Mouse CD276/B7-H3 mAb                                                                                                       |
| Remarks/Alias    | CD276; 4lg-B7-H3; B7-H3; B7H3; B7RP-2; CD276 antigen                                                                                          |
| Species          | Human                                                                                                                                         |
| GeneID (Human)   | 80381                                                                                                                                         |
| GeneID           | 80381                                                                                                                                         |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 29-245 of human CD276/B7-H3 (NP_079516.1). |
| Source           | Rabbit                                                                                                                                        |
| Category         | Monoclonal Antibodies                                                                                                                         |
| Application      | FC                                                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                                                  |
| SWISS            | Q5ZPR3-2                                                                                                                                      |
| Protein Weight   | 33kDa/52kDa/57kDa                                                                                                                             |
| Shipping         | Ice bag                                                                                                                                       |