

Abflo 647 Rabbit Anti-Human/Mouse Cd276/B7-H3 Monoclonal Antibody For Flow Cytometry - Q5Zpr3-2

Item Number: CM0011368

Introduction

Rabbit monoclonal antibody targeting CD276 (B7-H3) conjugated to ABflo 647, validated for flow cytometry on human and mouse specimens. Recombinant immunogen-based. UniProt Q5ZPR3-2. Ideal for immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse CD276/B7-H3 mAb
Remarks/Alias	B7H3; B7-H3; B7RP-2; 4lg-B7-H3
Species	Human
Gene ID (Human)	80381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-245 of human CD276 (NP_079516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q5ZPR3-2
Protein Weight	33kDa/52kDa/57kDa
Shipping	Ice bag