

Rabbit Anti-Human/Mouse Cd276/B7-H3 Monoclonal Antibody For Flow Cytometry - Q5Zpr3-2

Item Number: CM0008351

Introduction

-conjugated rabbit monoclonal antibody targeting CD276/B7-H3 (SWISS: Q5ZPR3-2). Suitable for flow cytometry with cross-reactivity to human and mouse. Ideal for immune checkpoint studies and tumor microenvironment analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse CD276/B7-H3 mAb
Remarks/Alias	CD276; 4lg-B7-H3; B7-H3; B7H3; B7RP-2; CD276 antigen
Species	Human
GeneID (Human)	80381
GeneID	80381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-245 of human CD276/B7-H3 (NP_079516.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q5ZPR3-2
Protein Weight	33kDa/52kDa/57kDa
Shipping	Ice bag