

## Rabbit Anti-Human/Mouse Cd276/B7-H3 Monoclonal Antibody For Flow Cytometry - Q5Zpr3-2

Item Number: CM0008003

### Introduction

-conjugated rabbit monoclonal antibody targeting human/mouse CD276/B7-H3, a costimulatory molecule. Validated for flow cytometry. Recognizes isoform 2 (Q5ZPR3-2). Suitable for immune checkpoint research.

[Learn More](#)

| Parameter        | Value                                                                                                                                        |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human/Mouse CD276/B7-H3 mAb                                                                                                      |
| Remarks/Alias    | B7H3; B7-H3; B7RP-2; 4Ig-B7-H3                                                                                                               |
| Species          | Human                                                                                                                                        |
| Gene ID (Human)  | 80381                                                                                                                                        |
| Gene ID          | 80381                                                                                                                                        |
| Immunogen        | Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 29-245 of human CD276/B7-H3 (NP_079516.1) |
| Source           | Rabbit                                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                                        |
| Application      | FC                                                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                                                 |
| SWISS            | Q5ZPR3-2                                                                                                                                     |
| Protein Weight   | 33kDa/52kDa/57kDa                                                                                                                            |
| Shipping         | Ice bag                                                                                                                                      |