

Rabbit Anti-Human/Mouse Cd351/Fcamr Mab - Q8Wwv6

Item Number: CM0008189

Introduction

High-affinity rabbit monoclonal antibody against CD351/FCAMR, conjugated to . Validated for flow cytometry with human and mouse samples. Ideal for studying IgA/IgM receptor-mediated immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse CD351/FCAMR mAb
Remarks/Alias	CD351; FCA/MR; FKSG87; Fcalpha/muR
Species	Human
GeneID (Human)	83953
GeneID	83953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-450 of human CD351/FCAMR (Q8WWV6-1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q8WWV6
Protein Weight	57kDa
Shipping	Ice bag