

Rabbit Anti-Human/Mouse Foxp3 Monoclonal Antibody For Fc (Intra) - Q9Bzs1

Item Number: CM0010070



Introduction

Rabbit anti-Human/Mouse FOXP3 monoclonal antibody for intracellular flow cytometry (FC intra). Cross-reacts with human and mouse FOXP3. Suitable for Treg research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse FOXP3 mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human/Mouse
GeneID (human)	50943
GeneID	50943
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human/mouse FOXP3 (NP_054728.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse
SWISS	Q9BZS1
Protein Weight	47 kDa/43 kDa/50 kDa/44 kDa
Shipping	Ice bag