

Rabbit Anti-Human/Mouse Integrin Alpha 5 (Itga5/Cd49E) Monoclonal Antibody For Flow Cytometry - P08648

Item Number: CM0004140



Introduction

-conjugated rabbit monoclonal antibody against human/mouse integrin alpha-5 (ITGA5, CD49e), validated for flow cytometry. Functions as fibronectin receptor involved in cell adhesion. Suitable for IVD research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse Integrin alpha 5 (ITGA5/CD49e) mAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A
Species	Human
Gene ID (Human)	3678
Gene ID	3678
Immunogen	Recombinant protein of human/mouse alpha 5 (ITGA5/CD49e)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag