

Rabbit Anti-Human/Mouse Pdgfra (Cd140A) Monoclonal Antibody For Fc - P16234

Item Number: CM0011547



Introduction

-conjugated rabbit monoclonal antibody targeting human/mouse PDGFR α (CD140a) extracellular domain. Validated for flow cytometry (FC). Reacts with human, mouse. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse PDGFR α /CD140a mAb
Remarks/Alias	PDGFRA; CD140A; PDGFR-2; PDGFR2; GAS9; RHEPDGFRA; platelet-derived growth factor receptor alpha
Species	Human
GeneID (human)	5156
GeneID	5156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human PDGFR α (NP_006197.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P16234
Protein Weight	24kDa/82kDa/122kDa
Shipping	Ice bag