

Rabbit Anti-Human/Mouse Pdgfra/Cd140A Monoclonal Antibody For Fc - P16234

Item Number: CM0010095



Introduction

-conjugated rabbit monoclonal antibody against human/mouse PDGFR α /CD140a, validated for flow cytometry. Suitable for studying PDGF signaling in development, wound healing, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse PDGFR α /CD140a mAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2
Species	Human
Gene ID (Human)	5156
Gene ID	5156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human PDGFR α (NP_006197.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P16234
Protein Weight	24kDa/82kDa/123kDa
Shipping	Ice bag