

Rabbit Anti-Human/Mouse/Monkey Foxp3 Monoclonal Antibody For Intracellular Flow Cytometry - Q9Bzs1

Item Number: CM0010523



Introduction

-conjugated rabbit monoclonal antibody specific for human, mouse, and cynomolgus FOXP3, validated for intracellular flow cytometry. Ideal for Treg research and immune homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse/Monkey FOXP3 mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human
GeneID (Human)	50943
GeneID	50943
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse, Cynomolgus
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag